

















Da:all

Products  
Information








# Da:all Products Line Up

## Da:all Premium Line

Brand Name	Model Number	Needle Type	Needle Gauge	Needle Length	USP	Thread (mm)	Folding	Shape
Da:all Premium Molding	CW18100-W3417	Blunt(W)	18G	100mm	3&4	170		
	CA18100-W3420	Blunt(A)	18G	100mm	3&4	200	folding	
	CA18100-W3417	Blunt(A)	18G	100mm	3&4	170		
	CA19100-W0120	Blunt(A)	19G	100mm	1	200	folding	
	CA19100-W0116	Blunt(A)	19G	100mm	1	160		
	ZA18100-W3420	Blunt(A)	18G	100mm	3&4	200	folding	
	ZA18100-W3417	Blunt(A)	18G	100mm	3&4	170		
	ZA19100-W0020	Blunt(A)	19G	100mm	0	200	folding	
ZA19100-W0016	Blunt(A)	19G	100mm	0	160			
Da:all Premium Emboss	KA18100-W3418	Blunt(A)	18G	100mm	3&4	180	folding	
	KA18100-W3416	Blunt(A)	18G	100mm	3&4	160		
	KA19100-W0118	Blunt(A)	19G	100mm	1	180	folding	
	KA19100-W0116	Blunt(A)	19G	100mm	1	160		
	KA18100-W3418	Blunt(L)	18G	100mm	3&4	180	folding	
KL19100-W0118	Blunt(L)	19G	100mm	1	180	folding		
Da:all Spiral Cobra	PAW19100-W0020	Blunt(W)	19G	100mm	0	200	folding	
	PAA19100-W0020	Blunt(A)	19G	100mm	0	200		
	PAL19100-W0020	Blunt(L)	19G	100mm	0	200		
	PAL2070-W2014	Blunt(L)	20G	70mm	2/0	140		
Da:all Spiral Double arm	AD11B00-W172	Blunt	21G	110mm	0	480	folding	
	AD11B00-W92	Blunt	21G	110mm	0	300		
	AD11B00-W72	Blunt	21G	80mm	0	240		
	AD16S00-W172	Spatula	21G	160mm	0	480		
	AD11S00-W172	Spatula	21G	110mm	0	480	folding	
	AD11S00-W92	Spatula	21G	110mm	0	300		
	AD08S00-W72	Spatula	21G	80mm	0	240		
	AD11S20-W172	Spatula	22G	110mm	2/0	480		
Da:all Braid	P8B1950L-W5070	Blunt(L)	19G	50mm	5/0 x 8	70	folding	
	P8B1938L-W5040	Blunt(L)	19G	38mm	5/0 x 8	40		
	P6B1950L-W7070	Blunt(L)	19G	50mm	7/0 x 16	70		
	P6B1938L-W7040	Blunt(L)	19G	38mm	7/0 x 16	40		
Da:all Premium Spring	PSS2160L-W3030	Blunt(L)	21G	60mm	3/0	30	folding	
	PSS2360L-W4030	Blunt(L)	23G	60mm	4/0	30		
	PSS2560L-W5030	Blunt(L)	25G	60mm	5/0	30		
	PSS2750L-W5020	Blunt(L)	27G	50mm	5/0	20		
	PSS2950L-W6020	Blunt(L)	29G	50mm	6/0	20		
Da:all Premium Braid soft	P2B2190L-W70100	Blunt(L)	21G	90mm	7/0 x 12	100	folding	
	P2B2160L-W7090	Blunt(L)	21G	60mm	7/0 x 12	90		
	P2B2160L-W7060	Blunt(L)	21G	60mm	7/0 x 12	60		
	P9B2390L-W70160	Blunt(L)	23G	90mm	7/0 x 9	160		
	P9B2390L-W70130	Blunt(L)	23G	90mm	7/0 x 9	130		
	P9B2360L-W7090	Blunt(L)	23G	60mm	7/0 x 9	90		
	P9B2360L-W7060	Blunt(L)	23G	60mm	7/0 x 9	60		

## Da:all Basic Line

Brand Name	Model Number	Needle Type	Needle Gauge	Needle Length	USP	Thread (mm)	Folding	Shape
Da:all Mono	M2590S-W40140	Sharp	25G	90mm	4/0	140	folding	
	M2760S-W5090	Sharp	27G	60mm	5/0	90		
	M2938S-W6050	Sharp	29G	38mm	6/0	50		
	M2925S-W6030	Sharp	29G	25mm	6/0	30		
Da:all Mono Screw	MS2590S-W40140	Sharp	25G	90mm	4/0	140	folding	
	MS2760S-W5090	Sharp	27G	60mm	5/0	90		
	MS2937S-W6050	Sharp	29G	38mm	6/0	50		
	MS2925S-W6030	Sharp	29G	25mm	6/0	30		
Da:all Double	D2660S-W7090	Sharp	26G	60mm	7/0 x 2	90	folding	
	D2638S-W7060	Sharp	26G	38mm	7/0 x 2	60		
Da:all Double Screw	DS2660S-W7090	Sharp	26G	60mm	7/0 x 2	90	folding	
	DS2638S-W7060	Sharp	26G	38mm	7/0 x 2	60		
Da:all I-Screw	MS3025U-W7030	Filler Needle	30G	25mm	7/0	30	folding	



## Leading technology company in advanced biotechnology, META-BIOMED



Established in 1990, META-BIOMED has been developing and producing raw materials for suture thread based on its distinguished technology accumulated over the last 28 years. In addition, META-BIOMED is specialized in biodegradable and non-absorbent thread and has acquired the own original technology of polymerizing the polymer required for various thread manufacturing.

- Over 70 intellectual property rights due to continuous research and development on biocompatible medical materials
- Developed differentiated global marketing with more than 250 sales networks in over 100 countries around the world
- Operating strict quality control system to ensure the highest quality product production
- Selected as 'World-class Products'(2003, 2007) and 'World-class 300 Companies' for innovative products(2015)
- Completed 'Meta Beauty R&D Center'(2018) to lead new technology development through bold R&D investment



## The first and only white PDO suture in Korea, Da:all



Please pay attention to the new brand value of Da:all with different dimension, differentiated technology, bold R&D investment and enthusiasm for the best.

- The only white PDO in Korea that distinguishes itself from other brands in violet color
- Transparent and undyed white PDO with no visibility during skin treatment
- Optimized for the thin layers of skin on the nose, forehead, under eyes of various races



# Worldwide First Undyed Lifting Thread



[Seoul Sales Office]

Tel +82-02-6933-7521~4 | Fax +82-02-6933-7525

E-mail : [cosmed@meta-biomed.com](mailto:cosmed@meta-biomed.com)

4F, 16, Gangnam-daero 95-gil, Seocho-gu, Seoul, Republic of Korea

<http://www.meta-biomed.com> | [www.facebook.com/DaallAesthetic](http://www.facebook.com/DaallAesthetic)

